

XDS-1 Inverted Biological Microscope



Microteb

Tel: 021 66 97 60 37

Fax: 021 66 90 33 54

www.microteb.com

Specification:

Order No.	Items	Standard Outfits
	Main Body (including left & right supporting arm plate, 4 piece of Allen screws & a wrench)	√
EW10	Eyepiece WF10×/ ϕ 20	√√
EW16	Eyepiece WF16×/ ϕ 14	√√
EC11	Centering Eyepiece DZ11×	√
EPr10b	Plan Article Eyepiece PF10×	√
B-OLPr10	LWD Plan Achromat Objective 10/0.25, 160/1.5, W.D. = 7.9mm	√
B-OLPr25	LWD Plan Achromat Objective 25/0.40, 160/1.5, W.D. = 5mm	√
B-OLPr40	LWD Plan Achromat Objective 40/0.60, 160/1.5, W.D. = 3mm	√
B-OLPnr10	LWD Plan Achromat Objective 10/0.25, 160/1.5, W.D. = 7.9mm	
B-OLPnr25	LWD Plan Achromat Objective 25/0.40, 160/1.5, W.D. = 5mm	√
B-OLPnr40	LWD Plan Achromat Objective 40/0.60, 160/1.5, W.D. = 3mm	
	35mm DF Camera With Shutter Control Line	
33X-4	Set for Stage Insert Plate	√
9J-8	Stage Insert plate (III)	√
9J-9	Stage Insert plate (IV)	√
33X-5	Stage Insert plate (V)	√
CL02	Stage Clips	√√
7X.5.1	Cultural Flask (bottom thickness 1.5mm)	√√
7X.5.2	Cultural Vessel (bottom thickness 1.5mm)	√√
RP06	Annular Spot 10×	
12X.2.3-19,20	Annular Spot 25×	√
RP07	Annular Spot 40×	
	Halogen Tungsten Lamp 6V20W (G4)	√√√
9J.7	Photo Adaptor With Photo Eyepiece S5×	√
33X. 3.2.1	Movable Condenser Insert Plate	√
TBR03	Seidentopf Binocular Head (inclined 45°)	√
	Micrometer (division 0.01mm)	
EPHp6.3	Photo Eyepiece S6.3×	
9J.15	Photo Eyepiece S0.65×For CCD	
CU43-3	C-mount	
	CCD TV & Record Attachment	